

Anti-CD66b/CEACAM8 antibody [G10F5] {APC} (STJA0002258)

STJA0002258

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Carcinoembryonic antigen-related cell adhesion molecule 8 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G10F5
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	1088
Gene Symbol	CEACAM8
Uniprot ID	CEAM8_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	