

## Anti-IgG1k antibody [MOPC-21] {PE/Cyanine5.5} (STJA0002340)

STJA0002340

### GENERAL INFORMATION

|                          |                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-IgG1, k is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                  |
| <b>Host/Source</b>       | Mouse                                                                                               |
| <b>Reactivity</b>        |                                                                                                     |

### PRODUCT PROPERTIES

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | MOPC-21                                                                                                      |
| <b>Concentration</b>       |                                                                                                              |
| <b>Conjugation</b>         | PE/Cyanine5.5                                                                                                |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA                                                                            |
| <b>Isotype</b>             | IgG1k                                                                                                        |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |
|---------------------------|
| <b>Gene ID</b>            |
| <b>Gene Symbol</b>        |
| <b>Uniprot ID</b>         |
| <b>Immunogen</b>          |
| <b>Immunogen Region</b>   |
| <b>Specificity</b>        |
| <b>Immunogen Sequence</b> |