

Anti-IgG2ak antibody [2A3] {PE/Cyanine7} (STJA0002453)
STJA0002453**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2a, k is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2A3
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence