

## Anti-IgG2bk antibody [LTF-2] {APC} (STJA0002475)

STJA0002475

### GENERAL INFORMATION

|                          |                                                                                                    |
|--------------------------|----------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                 |
| <b>Short Description</b> | Rat monoclonal antibody anti-IgG2b, k is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                 |
| <b>Host/Source</b>       | Rat                                                                                                |
| <b>Reactivity</b>        |                                                                                                    |

### PRODUCT PROPERTIES

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | LTF-2                                                                                                        |
| <b>Concentration</b>       | 0.2 mg/mL                                                                                                    |
| <b>Conjugation</b>         | APC                                                                                                          |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA                                                                            |
| <b>Isotype</b>             | IgG2bk                                                                                                       |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |
|---------------------------|
| <b>Gene ID</b>            |
| <b>Gene Symbol</b>        |
| <b>Uniprot ID</b>         |
| <b>Immunogen</b>          |
| <b>Immunogen Region</b>   |
| <b>Specificity</b>        |
| <b>Immunogen Sequence</b> |