

Anti-CD3/CD16/CD8a antibody [OKT3, 3G8, 5.1H11] {FITC/PE} (STJA0002509)
STJA0002509**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	NA monoclonal antibody anti-CD3-FITC/CD16+CD56-PE is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OKT3, 3G8, 5.1H11
Concentration	5 Mu L/Test
Conjugation	FITC/PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence