

## Anti-CD3/CD4/CD8a antibody [OKT3, RPA-T4, OKT8] {APC/FITC/PerCP-Cyanine5.5} (STJA0002510)

STJA0002510

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	NA monoclonal antibody anti-CD3-APC/CD4-FITC/CD8 Alpha-PerCP-Cyanine5.5 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	OKT3, RPA-T4, OKT8
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	APC/FITC/PerCP-Cyanine5.5
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>