

Anti-CD64 antibody [10.1] {EV450} (STJA0002530)

STJA0002530

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD64 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10.1
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA 2209
Gene Symbol	FCGR1BP FCGR1A
Uniprot ID	FCGRB_HUMAN FCGR1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	