

Anti-CD335/Ncr1 antibody [29A1.4] {PE/TR} (STJA0002561)

STJA0002561

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Natural cytotoxicity triggering receptor 1 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 29A1.4 |
| Concentration | 5 Mu L/Test |
| Conjugation | PE/TR |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------------------|-----------------------------|
| Gene ID | 17086 |
| Gene Symbol | Ncr1 |
| Uniprot ID | NCTR1_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |