

Anti-IgG2a, Lambda antibody [B39-4] {AF647} (STJA0002777)

STJA0002777

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2a, Lambda is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B39-4
Concentration	0.5 mg/mL
Conjugation	AF647
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2a, lambda
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence