

**Anti-Papillomavirus type 1 E2-protein antibody (208-218 aa) [1H10] (STJA0002812)**

STJA0002812

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Papillomavirus type 1 E2-protein (208-218 aa) is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
<b>Applications</b>	ELISA/WB/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	PV

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1H10
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein A purification
<b>Dilution Range</b>	ELISA 0.2 to 1 µg/ml; WB 0.2 to 1 µg/ml; IF 1 µg/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format.
<b>Formulation</b>	PBS with 50% glycerol
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	208-218 aa
<b>Immunogen</b>	208-218 aa
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	