

## Anti-Papillomavirus type 1 E2-protein antibody (310-410 aa) [3E8] (STJA0002813)

STJA0002813

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Papillomavirus type 1 E2-protein (310-410 aa) is suitable for use in ELISA and Immunofluorescence research applications.
<b>Applications</b>	ELISA/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	PV

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	3E8
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein A purification
<b>Dilution Range</b>	ELISA-0, 2-1 µg/ml; IF 0.33 to 20 µg/ml This is a conformational antibody and does not work in Western immunoblotting analysis or in other denaturative conditions. Monoclonal antibody working titer has to be established practically for each particula
<b>Formulation</b>	PBS with 50% glycerol
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	310-410 aa
<b>Immunogen Region</b>	310-410 aa
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081