

Anti-Papillomavirus type 1 E2-protein antibody (310-410 aa) [5H4] (STJA0002817)

STJA0002817

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Papillomavirus type 1 E2-protein (310-410 aa) is suitable for use in ELISA and Immunofluorescence research applications.
Applications	ELISA/IF
Host/Source	Mouse
Reactivity	PV

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5H4
Concentration	
Conjugation	Unconjugated
Purification	NA Protein A purification NA
Dilution Range	ELISA 0.1 to 0.2 µg/ml; IF 0.33 to 20 µg/ml This is a conformational antibody and does not work in Western immunoblotting analysis or in other denaturative conditions. Induces DNA/protein complex supershift in mobility shift assay.Monoclonal antibody
Formulation	PBS with 50% glycerol
Isotype	IgG
Storage	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	310-410 aa
Immunogen Region	310-410 aa
Specificity	
Immunogen Sequence	