

Anti-Hev b5 antibody [11] (STJA0002873)

STJA0002873

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Major latex allergen is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	HEVB5

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 164 000. Reacts as detection antibody in capture ELISA with Hev b5 binding monoclonal antibody A3-707-100
Formulation	PBS, 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Hev b5-MBP fusion protein expressed and purified from E. coli. Hev b5
Immunogen Region	
Specificity	Recombinant and native Hev b5
Immunogen Sequence	