

Anti-Hev b6.02 antibody [N3] (STJA0002875)

STJA0002875

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hevein is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Mouse
Reactivity	HEVB6, 02

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N3
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 132 000, WB 1. Reacts as detection antibody in capture ELISA with Hev b6.02 binding monoclonal antibody A3-709-100
Formulation	PBS, 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Hev b6.02-MBP fusion protein expressed and purified from E. coli. Hev b6.02
Immunogen Region	
Specificity	Recombinant and native Hev b6.02
Immunogen Sequence	