

Anti-DNase I antibody [4B10] (STJA0002876)

STJA0002876

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DNase I is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Mouse
Reactivity	DNASE1

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4B10
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.01 to 0.02 µg/ml; WB 0.1 to 0.2 µg/ml (bovine DNase I) Antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

