

## Anti-GFRa-4 antibody (STJA0002885)

STJA0002885

### GENERAL INFORMATION



<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Chicken polyclonal antibody anti-GFRa-4 is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
<b>Applications</b>	ELISA/WB/IF
<b>Host/Source</b>	Chicken
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	GFRa-4 affinity purified
<b>Dilution Range</b>	WB 0.1 to 0.2 µg/ml; ELISA 3 to 6.7 ng/ml; IF 0.3 to 1 µg/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
<b>Formulation</b>	PBS pH 7.4
<b>Isotype</b>	IgY
<b>Storage Instruction</b>	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	64096
<b>Gene Symbol</b>	GFRA4
<b>Uniprot ID</b>	GFRA4_HUMAN
<b>Immunogen</b>	Recombinant His-tagged human GFRa-4 protein produced using CHO-based proprietary suspension cell line. For production of hGFRa-4, glycosylphosphatidyl-inositol GPI-anchor was removed and protein was secreted to the cell culture supernatant. Protein w
<b>Immunogen Region</b>	
<b>Specificity</b>	Human GFRa-4, doesn't cross-react with hGFRa-1, hGFRa-2, hGFRa-3
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081