

Anti-S100-A8 protein antibody [5D8] (STJA0002913)

STJA0002913

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-S100-A8 protein is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5D8
Concentration	
Conjugation	Unconjugated
Purification	NA Protein G purification NA
Dilution Range	ELISA 2 µg/ml
Formulation	PBS with 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6279
Gene Symbol	S100A8
Uniprot ID	S10A8_HUMAN
Immunogen	Human S100-A8
Immunogen Region	
Specificity	Human S100-A8
Immunogen Sequence	