

Anti-Recombinant-MANF antibody [3D4] (STJA0002919)

STJA0002919

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-MANF is suitable for use in ELISA, Immunofluorescence and Immunohistochemistry-C research applications.
Applications	ELISA/IF/IHC-C
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3D4
Concentration	
Conjugation	Unconjugated
Purification	MabSelect affinity chromatography
Dilution Range	ELISA 0.02 to 0.1 µg/ml; IF 1 to 20 µg/ml; IHC-C 5 µg/ml. Conformational antibody, not suitable for Western Blot application. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG2a
Storage	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
Instruction	

TARGET INFORMATION

Gene ID	7873
Gene Symbol	MANF
Uniprot ID	MANF_HUMAN
Immunogen	Recombinant human MANF protein produced using CHO-based proprietary suspension cell line. Immunogen is purified from cell culture supernatant. Human MANF
Immunogen Region	
Specificity	Human MANF
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081