

Anti-Recombinant-SARS-CoV-2 Nucleocapsid protein antibody [82C3] (STJA0002929)

STJA0002929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Cho Cells monoclonal antibody anti-Recombinant-SARS-CoV-2 Nucleocapsid protein is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	CHO cells
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	82C3
Concentration	
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Dilution Range	Indirect ELISA 0.005 – 0.01 ng/ml
Formulation	PBS pH 7.4
Isotype	IgG1 kappa
Storage Instruction	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein produced by CHO-based proprietary suspension cell line.
Immunogen Region	
Specificity	SARS-CoV-2 Nucleocapsid protein
Immunogen Sequence	