

Anti-Recombinant-interleukin-6 antibody [11C5] (STJA0002936)

STJA0002936

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Cho Cells monoclonal antibody anti-Recombinant-interleukin-6 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	CHO cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	11C5
Concentration	
Conjugation	Unconjugated
Purification	NA Protein A affinity chromatography following desalting NA
Dilution Range	
Formulation	PBS pH 7.4
Isotype	IgG1
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
Uniprot ID	IL6_HUMAN
Immunogen	Recombinant hIL-6 protein produced by CHO-based cell line (expressed by QMCF Technology).
Immunogen Region	
Specificity	Reacts with human IL-6.No cross-reactivity with IL-1 Alpha, IL-1 Beta, IL-2, IL-3, IL-4, IL-8, IFN Gamma and TNF Alpha proteins.
Immunogen Sequence	