

**Anti-Recombinant-Cardiac Troponin I antibody [5F8] (STJA0002946)**

STJA0002946

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Cho Cells monoclonal antibody anti-Recombinant-Cardiac Troponin I is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	CHO cells
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5F8
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	MabSelect affinity chromatography following gel filtration
<b>Dilution Range</b>	
<b>Formulation</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> (SPTN)
<b>Isotype</b>	IgG2 lambda
<b>Storage Instruction</b>	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Recombinant human cardiac troponin complex protein produced by CHO-based proprietary suspension cell line.
<b>Immunogen Region</b>	
<b>Specificity</b>	human cardiac troponin I
<b>Immunogen Sequence</b>	