

Anti-Recombinant-HEL Mouse IgG2a antibody (STJA0002952)

STJA0002952

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Cho Cells monoclonal antibody anti-Recombinant-HEL Mouse IgG2a is suitable for use in research applications.
Applications	
Host/Source	CHO cells
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Hitrap Mabselect Sure + Superdex 200 Increase
Dilution Range	
Formulation	PBS pH 7.4
Isotype	
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	Lysozyme from chicken egg white
Immunogen Sequence	