

Anti-Chrnb4 antibody (Cytoplasmic Loop) (STJA0003704)

STJA0003704

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-nAChR b4 (Cytoplasmic Loop) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000
Formulation	100 µl in 10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 µg per ml BSA and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	108015
Gene Symbol	Chrnb4
Uniprot ID	ACHB4_MOUSE
Immunogen	Fusion protein from the cytoplasmic loop of the beta 4 subunit of rat nAChR.
Immunogen Region	Cytoplasmic Loop
Specificity	
Immunogen Sequence	



Western blot of mouse habenula lysate showing specific immunolabeling of the ~52 kDa nAChR Beta 4 protein.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081