

Anti-NMDAR1 C2 antibody (STJA0003708)

STJA0003708

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NMDAR1 C2 is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Rabbit
Reactivity	Mouse/Rat/Bovine/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from pooled serum.
Dilution Range	WB 1:1000 IHC 1:1000-1:2000
Formulation	Lyophilized
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24408
Gene Symbol	Grin1
Uniprot ID	NMDZ1_RAT
Immunogen	Synthetic peptide corresponding to amino acid residues specific to the NR1 subunit, C2 splice variant insert conjugated to KLH.
Immunogen Region	
Specificity	
Immunogen Sequence	