

Anti-Nuclei antibody [235-1] (STJA0003728)

STJA0003728

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Nuclei is suitable for use in Immunohistochemistry and Immunocytochemistry research applications.
Applications	IHC/ICC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	235-1
Concentration	
Conjugation	Unconjugated
Purification	This antibody was concentrated tissue culture from supernatant.
Dilution Range	IHC 1:100
	ICC 1:100
Formulation	100 ul
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human cell homogenate.
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunostaining of HeLa cells showing specific labeling of nuclei using the anti-nuclei antibody (STJA0003728, 1:100, green).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081