

Anti-POSTN-Pan antibody (Fascilin Domain) (STJA0003746)

STJA0003746

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Periostin-Pan (Fascilin Domain) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Bovine/Canine/Feline/Finch/Fish/Goat/Guinea Pig/Hamster/Horse/Non-Human Primates/Rabbit/Sheep/Vole/Xenopus/Zebrafish

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	This antibody was antigen affinity purified from serum.
Dilution Range	WB 1:1000 IHC 1:100
Formulation	100 ul in PBS
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	50706
Gene Symbol	Postn
Uniprot ID	POSTN_MOUSE
Immunogen	Synthetic peptide corresponding to amino acid residues from the fascilin domain 1 of mouse periostin.
Immunogen Region	Fascilin Domain
Specificity	
Immunogen Sequence	



Western blot of rat lung lysate showing specific immunolabeling of the ~93 kDa periostin protein triplet.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081