

Anti-Vimentin antibody [LN6] (STJA0003872)

STJA0003872

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Vimentin is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LN6
Concentration	
Conjugation	Unconjugated
Purification	Protein L purified culture supernatant
Dilution Range	WB 1:1000 IHC 1:100
Formulation	100 ul in PBS
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
Uniprot ID	VIME_HUMAN
Immunogen	Cell nuclei preparation from adult human thymus.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot of HeLa cell lysate showing specific immunolabeling of the ~50 kDa vimentin protein.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081