

Anti-CD3 antibody [17A2] {EV450} (STJA0003893)

STJA0003893

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	0.5 mg/mL
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12501 12503 12500 Cd3e Cd247 CD3E_MOUSE CD3Z_MOUSE CD3D_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	