

Anti-CD11b antibody [M1/70] {ER780} (STJA0003901)

STJA0003901

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	0.5 mg/mL
Conjugation	ER780
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684
	NA
Gene Symbol	ITGAM
	Itgam
Uniprot ID	ITAM_HUMAN
	ITAM_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	