

Anti-Cd4 antibody [GK1.5] {EV450} (STJA0003904)

STJA0003904

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GK1.5
Concentration	0.5 mg/mL
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
Uniprot ID	CD4_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	