

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] {ER780} (STJA0003919)
STJA0003919**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD45R/B220 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RA3.3A 1/6.1
Concentration	0.5 mg/mL
Conjugation	ER780
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	5788 19264
Gene Symbol	PTPRC Ptprc
Uniprot ID	PTPRC_HUMAN PTPRC_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	