

## Anti-Ly6C antibody [Monts 1] {EV450} (STJA0003924)

STJA0003924

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Ly6C is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	Monts 1
<b>Concentration</b>	5 µg L/Test
<b>Conjugation</b>	EV450
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">17067</a>
	<a href="#">100041546</a>
<b>Gene Symbol</b>	<a href="#">Ly6c1</a>
	<a href="#">Ly6c2</a>
<b>Uniprot ID</b>	<a href="#">LY6C1_MOUSE</a>
	<a href="#">LY6C2_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	