

Anti-Ly6C antibody [Monts 1] {ER780} (STJA0003927)

STJA0003927

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly6C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Monts 1
Concentration	0.5 mg/mL
Conjugation	ER780
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	100041546
	17067
Gene Symbol	Ly6c2
	Ly6c1
Uniprot ID	LY6C2_MOUSE
	LY6C1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	