

**Anti-CD40 antibody [G28.5] {ER780} (STJA0003934)**

STJA0003934

**GENERAL INFORMATION**

|                          |                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                              |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 5 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                                              |
| <b>Host/Source</b>       | Mouse                                                                                                                                           |
| <b>Reactivity</b>        | Human                                                                                                                                           |

**PRODUCT PROPERTIES**

|                       |                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Monoclonal                                                                                                   |
| <b>Clone ID</b>       | G28.5                                                                                                        |
| <b>Concentration</b>  | 5 Mu L/Test                                                                                                  |
| <b>Conjugation</b>    | ER780                                                                                                        |
| <b>Purification</b>   |                                                                                                              |
| <b>Dilution Range</b> |                                                                                                              |
| <b>Formulation</b>    | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>        | IgG1k                                                                                                        |
| <b>Storage</b>        | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |
| <b>Instruction</b>    |                                                                                                              |

**TARGET INFORMATION**

|                    |                            |
|--------------------|----------------------------|
| <b>Gene ID</b>     | <a href="#">958</a>        |
| <b>Gene Symbol</b> | <a href="#">CD40</a>       |
| <b>Uniprot ID</b>  | <a href="#">TNR5_HUMAN</a> |
| <b>Immunogen</b>   |                            |
| <b>Immunogen</b>   |                            |
| <b>Region</b>      |                            |
| <b>Specificity</b> |                            |
| <b>Immunogen</b>   |                            |
| <b>Sequence</b>    |                            |