

Anti-CD34 antibody [RAM34] {APC} (STJA0003941)

STJA0003941

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD34 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RAM34
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12490
Gene Symbol	Cd34
Uniprot ID	CD34_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	