

Anti-Cd4 antibody (Domain 1) [OX-38] {EV450} (STJA0003944)

STJA0003944

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD4 (Domain 1) is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-38
Concentration	0.5 mg/mL
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24932
Gene Symbol	Cd4
Uniprot ID	CD4_RAT
Immunogen	Domain 1
Immunogen Region	Domain 1
Specificity	
Immunogen	
Sequence	