

Anti-CD3 antibody [G4.18] {PE/TR} (STJA0003967)

STJA0003967

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G4.18
Concentration	5 µg L/Test
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG3k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	25710 300678
Gene Symbol	Cd3d Cd3g
Uniprot ID	CD3D_RAT CD3G_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	