

**Anti-CD3 antibody [G4.18] {PE/TR} (STJA0003968)**

STJA0003968

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	G4.18
<b>Concentration</b>	0.2 mg/mL
<b>Conjugation</b>	PE/TR
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG3k
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">25710</a> <a href="#">300678</a>
<b>Gene Symbol</b>	<a href="#">Cd3d</a> <a href="#">Cd3g</a>
<b>Uniprot ID</b>	<a href="#">CD3D_RAT</a> <a href="#">CD3G_RAT</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	