

**Anti-CD56/NCAM1 antibody [5.1H11] {FITC} (STJA0003976)**

STJA0003976

**GENERAL INFORMATION**

|                          |                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                          |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Neural cell adhesion molecule 1 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                          |
| <b>Host/Source</b>       | Mouse                                                                                                                       |
| <b>Reactivity</b>        | Human                                                                                                                       |

**PRODUCT PROPERTIES**

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | 5.1H11                                                                                                       |
| <b>Concentration</b>       | 5 Mu L/Test                                                                                                  |
| <b>Conjugation</b>         | FITC                                                                                                         |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>             | IgG1k                                                                                                        |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|                         |                             |
|-------------------------|-----------------------------|
| <b>Gene ID</b>          | <a href="#">4684</a>        |
| <b>Gene Symbol</b>      | <a href="#">NCAM1</a>       |
| <b>Uniprot ID</b>       | <a href="#">NCAM1_HUMAN</a> |
| <b>Immunogen</b>        |                             |
| <b>Immunogen Region</b> |                             |
| <b>Specificity</b>      |                             |
| <b>Immunogen</b>        |                             |
| <b>Sequence</b>         |                             |