

Anti-CD56/NCAM1 antibody [5.1H11] {EV450} (STJA0003982)

STJA0003982

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neural cell adhesion molecule 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5.1H11
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	4684
Gene Symbol	NCAM1
Uniprot ID	NCAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	