

Anti-CD117/Kit antibody [2B8] {PE/Cyanine7} (STJA0003986)

STJA0003986

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Mast/stem cell growth factor receptor Kit is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2B8
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16590
Gene Symbol	Kit
Uniprot ID	KIT_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	