

Anti-CD27 antibody [O323] {ER780} (STJA0003997)

STJA0003997

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD27 antigen is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | O323 |
| Concentration | 5 Mu L/Test |
| Conjugation | ER780 |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 939 |
| Gene Symbol | CD27 |
| Uniprot ID | CD27_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |