

Anti-CD45/Ptpnc antibody [30-F11] {PE/TR} (STJA0004000)

STJA0004000

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	30-F11
Concentration	5 Mu L/Test
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	19264
Gene Symbol	Ptpnc
Uniprot ID	PTPRC_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	