

Anti-CD45.1/Ptprc antibody [A20] {EV450} (STJA0004007)

STJA0004007

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | A20 |
| Concentration | 0.5 mg/mL |
| Conjugation | EV450 |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------------------|-----------------------------|
| Gene ID | 19264 |
| Gene Symbol | Ptprc |
| Uniprot ID | PTPRC_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |