

## Anti-PDCD1 antibody [EH12.2H7] {PerCP/Cyanine5.5} (STJA0004013)

STJA0004013

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Programmed cell death protein 1 is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	EH12.2H7
<b>Concentration</b>	5 Mu L/Test
<b>Conjugation</b>	PerCP/Cyanine5.5
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	5133
<b>Gene Symbol</b>	PDCD1
<b>Uniprot ID</b>	PDCD1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	