

Anti-CD57 antibody [HNK-1] {EV450} (STJA0004020)
STJA0004020**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD57 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HNK-1
Concentration	5 Mu L/Test
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	27087
Gene Symbol	B3GAT1
Uniprot ID	B3GA1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	