

Anti-HLA-DR antibody [L243] {EV450} (STJA0004025)

STJA0004025

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-DR is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	L243
Concentration	5 µg L/Test
Conjugation	EV450
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3122 3123
Gene Symbol	HLA-DRA HLA-DRB1
Uniprot ID	DRA_HUMAN DRB1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	