

Anti-CD326 antibody [G8.8] {AF488} (STJA0004026)

STJA0004026

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Epithelial Cell Adhesion Molecule is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G8.8
Concentration	5 Mu L/Test
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17075
Gene Symbol	Epcam
Uniprot ID	EPCAM_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	