

Anti-CD64 antibody [X54-5/7.1] {ER780} (STJA0004032)

STJA0004032

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-High Affinity Immunoglobulin Gamma Fc Receptor I is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	X54-5/7.1
Concentration	5 Mu L/Test
Conjugation	ER780
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	14129
Gene Symbol	Fcgr1
Uniprot ID	FCGR1_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	