

Anti-IL-17A antibody [17F3] {PE/Cyanine7} (STJA0004062)

STJA0004062

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-17a is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17F3
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16171
Gene Symbol	IL17a
Uniprot ID	IL17_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	